

# GC Ortho **Connect**<sup>TM</sup>



Light Cure Orthodontic Adhesive  
**One Step System**

**GC Orthodontics**  
GLOBAL EXPERTS IN ORTHODONTICS


# GC Ortho **Connect™** 90-1100-0000 (2 syringes of 2,9 g each)

## Light Cure Orthodontic Adhesive


After 18 months of successful research trials in Japan, GC Corporation is proud to launch GC Ortho Connect\*, a new orthodontic light cure adhesive for brackets and bondable buccal tubes. Now available for the European orthodontic market through GC Orthodontics.



When using GC Ortho Connect, **no primer is necessary**, since it is already incorporated into the adhesive, thus creating a **simplified bonding process** that **reduces costs**.



The **adhesive's ideal viscosity** insures reliable and precise bracket placement.




An ergonomic syringe guarantees **accurate dosage**, providing a practical and economic adhesive solution. The syringe's cap prevents dry out and thus decreasing potential waste.




**Excess adhesive is detected by a black light**, which aids in fast and efficient cleanup.

## Highly Resistant



### Excellent Bond Strength

Featuring superior flexural strength - fundamental in supporting mastication pressure.



### Reduced Bond Failure

The low elasticity component of GC Ortho Connect, compared to other composites in the industry, allows for high impact resistance without bond failure.

GC Ortho Connect provides reliable and versatile bonding on metal braqué and buccal tubes, ceramic and plastic brackets.



## GC Ortho **Etching Gel** 90-1400-0000

2 syringes (4,8 g each)  
40 dispensing tips

GC Orthodontics Europe GmbH  
Harkortstrasse 2  
D-58339 Breckerfeld  
Germany

