

# Technology that'll make a splash in your practice.



- **Eliminates water contaminants**  
Vaporizes tap water to kill and remove 99% of most contaminants.
- **Protects your equipment**  
Aquastat produces the highly pure water necessary for autoclave sterilization in your office.
- **Fast distilling/energy efficient**  
Aquastat processes 3.78 litres (1 gallon) of distilled water in just 4.25 hours with very low power consumption.
- **Safe, clean water from any source**  
Whether your water source is spring fed, municipally piped, soft or hard, the Aquastat assures you of fresh, safe, clean water every time.
- **Portable and easy-to-use**  
No special wiring or plumbing required – just plug it in! And Aquastat is easily portable at just 4.5 kg (10 lb.).
- **Removable boiling chamber**  
For ease of use and cleaning.
- **Special heating base**  
A uniquely designed heating base reduces scale buildup.
- **New design**  
A new contemporary design compliments any office decor.

**AQUASTAT**

# AQUASTAT

The cleanest water in a sea of technology



- Vaporization
- Condensation
- 99% Free of most contaminants

Weight: 4.5 kg (10 lb.)

Dimensions: 40.6 x 24.5 x 38.4 cm (L x W x H)  
(16 x 9 5/8 x 15 1/8")

## A new wave breakthrough provides clean, fresh water

Protect your equipment by using only steam-process distilled water in your practice. SciCan's AQUASTAT distiller allows you to forget about buying expensive bottled water or risking damage to your equipment from poor quality distilled water. A patented energy saving design with unique distillation technology kills harmful bacteria and viruses by vaporizing water and then reconstituting it through condensation. The result is water that's 99% free of most contaminants.

AQUASTAT processes soft or hard water from any source including those that are spring fed or municipally piped. Simply plug it in, fill with tap water and switch on.

With its energy efficient technology, AQUASTAT processes 3.78 litres (1 gallon) of water in just 4.25 hours. Combine this fast process with a compact, light-weight design that doesn't require any special plumbing or wiring – and you'll have all the fresh water you need.



### SciCan Ltd

1440 Don Mills Road, Toronto, Ontario M3B 3P9 CANADA  
Phone (416) 445-1600; Fax (416) 445-2727;  
Toll-free 1-800-667-7733

### SciCan GmbH

Kurzes Geländ 10 D-86156 Augsburg GERMANY  
Phone: +49 (0) 821 56 74 56-0; Fax: +49 (0) 821 56 74 56-99

### SciCan Inc

701 Technology Drive, Canonsburg, PA 15317 USA  
Phone (724) 820-1600; Fax (724) 820-1479;  
Toll-free 1-800-572-1211

### SciCan Medtech

Alpenstrasse 16 6300 Zug SWITZERLAND  
Local: (41- 41) 727.70.27; Fax: (41- 41) 727.70.29

[www.scican.com](http://www.scican.com)

**SciCan** Dental

Your Infection Control Specialist™